





MAINTENANCE TRACKING TOOL

Pre-Maintenance

Usage/Calendar based maintenance (UBM/CBM) NOTE! Refer to the PETtrace Service Manual - Maintenance (direction 2169049-100) for detailed instructions, apply LOTO and use PPE.	
System ID:	NE209962
Maintenance performed in accordance with instructions as outlined in the PETtrace Service Manual - Maintenance (direction 2169049-100) (signature (typed and signed)): 	

Location	Action	Labor time (min.)	Sign.	For only optional operations note down if the operation is performed or not															
Vacuum	NOTE! Hydrogen gas flow should be on as for normal production. <ul style="list-style-type: none"> • Read and record the vacuum pressure • Perform a BEV leak check : open the BEV for 2 minutes and close it. After 10 minutes open again the BEV,the vacuum value must not reach the value of 1.0*E-5 																		
	<table border="1"> <tr> <td>Vacuum pressure readout</td> <td>Gas flow(sccm):</td> <td>4.0</td> </tr> </table>				Vacuum pressure readout	Gas flow(sccm):	4.0												
	Vacuum pressure readout				Gas flow(sccm):	4.0													
	<table border="1"> <thead> <tr> <th>Gauge number</th> <th>Pressure (x10-) without gas</th> <th>Pressure (x10-) with gas</th> </tr> </thead> <tbody> <tr> <td>A1 (4 on TCS 1001):</td> <td>3,60E-03</td> <td>3,50E-02</td> </tr> <tr> <td>A2 (13 on TCS 1001):</td> <td>UR</td> <td>UR</td> </tr> <tr> <td>B1 (14 on TCS 1001):</td> <td>5,30E-08</td> <td>1,20E-05</td> </tr> </tbody> </table>				Gauge number	Pressure (x10-) without gas	Pressure (x10-) with gas	A1 (4 on TCS 1001):	3,60E-03	3,50E-02	A2 (13 on TCS 1001):	UR	UR	B1 (14 on TCS 1001):	5,30E-08	1,20E-05			
	Gauge number				Pressure (x10-) without gas	Pressure (x10-) with gas													
	A1 (4 on TCS 1001):				3,60E-03	3,50E-02													
	A2 (13 on TCS 1001):				UR	UR													
	B1 (14 on TCS 1001):				5,30E-08	1,20E-05													
	<table border="1"> <thead> <tr> <th colspan="3">TPG parameters</th> </tr> <tr> <th></th> <th>Low limit (x10-)</th> <th>High limit (x10-)</th> </tr> </thead> <tbody> <tr> <td>A1:</td> <td>1,00E-01</td> <td>7,00E-02</td> </tr> <tr> <td>A2:</td> <td>1,80E-05</td> <td>7,00E-01</td> </tr> <tr> <td>B1:</td> <td>2,00E-01</td> <td>3,00E-05</td> </tr> </tbody> </table>				TPG parameters				Low limit (x10-)	High limit (x10-)	A1:	1,00E-01	7,00E-02	A2:	1,80E-05	7,00E-01	B1:	2,00E-01	3,00E-05
	TPG parameters																		
	Low limit (x10-)	High limit (x10-)																	
A1:	1,00E-01	7,00E-02																	
A2:	1,80E-05	7,00E-01																	
B1:	2,00E-01	3,00E-05																	

Vacuum	• Press OFF on the VCU, followed by VENT, read and record the current VENT time			
	Vacuum VENT time			
	System software			
	Subsystem	Version		
	Master:	3,6		
	ACS:	4.3.2		
	Service System:	3.6.0		
	Manager:	TSA		
	Informix (only applicable to SUN-Master Station):	NA		
	VENT time:	2022-11-06 16:13		

Comments:	comments
<p>PHOTO:</p>	<p style="text-align: center;">Photo name:</p> 
	<p style="text-align: center;">Photo name:</p> 
	<p style="text-align: center;">Photo name:</p> 