





**MAINTENANCE TRACKING TOOL**

**Pre-Maintenance**

<b>Usage/Calendar based maintenance (UBM/CBM)</b> <b>NOTE! Refer to the PETtrace Service Manual - Maintenance (direction 2169049-100) for detailed instructions, apply LOTO and use PPE.</b>	
System ID:	NE209962
Maintenance performed in accordance with instructions as outlined in the PETtrace Service Manual - Maintenance (direction 2169049-100) (signature (typed and signed)): 	

Location	Action	Labor time (min.)	Sign.	For only optional operations note down if the operation is performed or not															
Vacuum	<b>NOTE! Hydrogen gas flow should be on as for normal production.</b> <ul style="list-style-type: none"> <li>• Read and record the vacuum pressure</li> <li>• Perform a BEV leak check : open the BEV for 2 minutes and close it. After 10 minutes open again the BEV,the vacuum value must not reach the value of 1.0*E-5</li> </ul>																		
	<table border="1"> <tr> <td><b>Vacuum pressure readout</b></td> <td><b>Gas flow(sccm):</b></td> <td>4.0</td> </tr> </table>				<b>Vacuum pressure readout</b>	<b>Gas flow(sccm):</b>	4.0												
	<b>Vacuum pressure readout</b>				<b>Gas flow(sccm):</b>	4.0													
	<table border="1"> <thead> <tr> <th>Gauge number</th> <th>Pressure (x10-) without gas</th> <th>Pressure (x10-) with gas</th> </tr> </thead> <tbody> <tr> <td>A1 (4 on TCS 1001):</td> <td>3,60E-03</td> <td>3,50E-02</td> </tr> <tr> <td>A2 (13 on TCS 1001):</td> <td>UR</td> <td>UR</td> </tr> <tr> <td>B1 (14 on TCS 1001):</td> <td>5,30E-08</td> <td>1,20E-05</td> </tr> </tbody> </table>				Gauge number	Pressure (x10-) without gas	Pressure (x10-) with gas	A1 (4 on TCS 1001):	3,60E-03	3,50E-02	A2 (13 on TCS 1001):	UR	UR	B1 (14 on TCS 1001):	5,30E-08	1,20E-05			
	Gauge number				Pressure (x10-) without gas	Pressure (x10-) with gas													
	A1 (4 on TCS 1001):				3,60E-03	3,50E-02													
	A2 (13 on TCS 1001):				UR	UR													
	B1 (14 on TCS 1001):				5,30E-08	1,20E-05													
	<table border="1"> <thead> <tr> <th colspan="3">TPG parameters</th> </tr> <tr> <th></th> <th>Low limit (x10-)</th> <th>High limit (x10-)</th> </tr> </thead> <tbody> <tr> <td>A1:</td> <td>1,00E-01</td> <td>7,00E-02</td> </tr> <tr> <td>A2:</td> <td>1,80E-05</td> <td>7,00E-01</td> </tr> <tr> <td>B1:</td> <td>2,00E-01</td> <td>3,00E-05</td> </tr> </tbody> </table>				TPG parameters				Low limit (x10-)	High limit (x10-)	A1:	1,00E-01	7,00E-02	A2:	1,80E-05	7,00E-01	B1:	2,00E-01	3,00E-05
	TPG parameters																		
	Low limit (x10-)	High limit (x10-)																	
A1:	1,00E-01	7,00E-02																	
A2:	1,80E-05	7,00E-01																	
B1:	2,00E-01	3,00E-05																	

Vacuum	• Press OFF on the VCU, followed by VENT, read and record the current VENT time			
	Vacuum VENT time			
	System software			
	Subsystem	Version		
	Master:	3,6		
	ACS:	4.3.2		
	Service System:	3.6.0		
	Manager:	TSA		
	Informix (only applicable to SUN-Master Station):	NA		
	VENT time:	2022-11-06 16:13		

Comments:	comments
<p>PHOTO:</p>	<p style="text-align: center;"><b>Photo name:</b></p> 
	<p style="text-align: center;"><b>Photo name:</b></p> 
	<p style="text-align: center;"><b>Photo name:</b></p> 