NIGHT INCIDENT



Report

		A Novartis Company
Site	Barcelona	
Date and time	2024-01-22 00:38	
Attended		
Phone Site Site		
Department		
Engineering		
Scope		
Yes No No		
System		
Distribution - Theodorico		
System description		
Describe the incident		
0:38 When opening the autoclave chamber, water comes ou Technicians are asked to clean the water and test the Theo a I ask you to tilt the sensor with a flat screwdriver so that it b	nd auto05 1:09 Theo does not read the	e vial after filling.
Solved?		
Yes No No		
Derived?		
Yes No X		
To whom?		
Yes - Site Yes - Production manager manager	Yes - QA Yes - QC QI	()thers
Who is other?		
Production affected?		
No X Yes - Totally cancelled Yes - Partially can	ncelled Yes - Delayed	
Any comment?		
If possible, describe the solution		
Cleaning the Theo hot cell and tilting the vial sensor		
If possible, describe the root cause		
drainage problems in the theodorico autoclave and crashes i	n the robot	
TotalTimeOfIntervention		