NIGHT INCIDENT Report	Advanced Accelerator Applications ANovantic Company
Site	Oporto
Date and time	2024-08-22 02:06
Attended Phone Site 🗙	
Department	
Engineering	
Scope	
Yes X No	
System	
Synthesizers - TRASIS AIO	
System description	

Describe the incident

Error in step 24 (FDOPA synthesis) with 2 different cassettes - error regarding vacuum on the water recovery vial. After performing tests to discard the module, it was used a 3rd cassette and no error again. Consumable problem.

Solved?
Yes X No
Derived?
Yes No 🗙
To whom?
Yes - Site managerYes - Production managerYes - QAYes - QCYes - QPOthers
Who is other?
Production affected?
No
Any comment?
If possible, describe the solution
If possible, describe the root cause
TotalTimeOfIntervention
60.0